

# Identification of Genes Induced by Taxol Application Using a Combination of Differential Display RT-PCR and DNA Microarray Analysis

Ei-ichiro Fukusaki<sup>a</sup>, \* Takashi Oishi<sup>a</sup>, Hozumi Tanaka<sup>b</sup>, Shin-ichiro Kajiyama<sup>a</sup> and Akio Kobayashi<sup>a</sup>

<sup>a</sup> Department of Biotechnology, Graduate School of Engineering, Osaka University,  
2-1 Yamadaoka, Suita, 565-0871, Japan  
Fax: +81-6-6879-7424. E-mail: fukusaki@bio.eng.osaka-u.ac.jp

<sup>b</sup> Medical Products Sector, Nitto Denko Corporation, 1-1-2 Shimohozumi, Ibaraki, 567-8680, Japan

\*Author for correspondance and reprint requests

Z. Naturforsch. **56c**, 814-819 (2001); received March 26/May 22, 2001

Taxol, HeLa Cells, DNA Microarray

The differential display reverse transcriptional polymerase chain reaction (DD-RT-PCR) was used to hunt for cDNA fragments specifically expressed by taxol treatment of HeLa cells. Forty-eight cDNA clones were differentially displayed through the experiments. The cDNA fragments obtained were separately spotted onto glass slides to prepare a tailor-made DNA chip. The gene expression pattern of differentially displayed cDNA fragments were checked by DNA microarray analysis.